



Specialised  
Orthotic  
Services

**POZIFORM®**

**GROIN HARNESS**

*Installation &  
User Instructions*

**POZIFORM®** is a registered trademark of Specialised Orthotic Services Ltd.

Thank you for purchasing your **POZIFORM® GROIN HARNESS**.

These harnesses are available in a range of sizes to suit most users everyday needs, designed to help fix the position of the pelvis without putting any pressure on the abdominal area.

All **POZIFORM®** products are supplied with a comprehensive 3 months warranty from the date of purchase (For terms and conditions please consult our Customer Services Department).

Please note that the groin harness should be used in conjunction with the **FITTINGS** supplied and **must not** be used as a transportation safety belt in a moving vehicle.

To fit your **POZIFORM®** Groin Harness please follow these instructions supplied.

## Specialised Orthotic Services Ltd

Units 127/128 Fauld Ind Park

Tutbury, Nr Burton-on-Trent,

Staffordshire. ENGLAND. DE13 9HS

Tel: 0044 (0)1283 520 400 Fax: 0044 (0)1283 520 401

E-mail: [enquiries@specialorthotic.com](mailto:enquiries@specialorthotic.com)


Web: [specialisedorthoticservices.co.uk](http://specialisedorthoticservices.co.uk)



Standard Groin Harness

### KIT LIST & TOOL REQUIREMENTS

**Your Harness Kit Consists Of The Following:**



1 x Fitting Instructions  
1 x Groin Harness

**Tools Required (not supplied):**

NO TOOLS REQUIRED

## FITTING INSTRUCTIONS

**BEFORE** fitting the harness to any wheelchair/seating system please consult the manufacturers hand book (or contact them direct) to find out if it is safe to do so.

### 1) Fitting to Wheelchair or Special Seating



- 1.1) a) Place the groin harness on the seat base with the straps on the underside.  
b) Wrap the two straps (with the male & female buckles) around the back of the wheelchair uprights and clip together. Pull the two straps to tension.

### 2) User Guide: The following user instructions are applicable when fitted to the wheelchair or Seating System.



- 2.1) Once the user is sat on the groin harness in the wheelchair wrap the harness pads around the thighs (as shown), Grasp the two halves of the buckle and securely press together. Finally tension the Harness by grasping the buckle and pulling the webbing through the buckle until the harness is not too loose to be ineffective nor too tight to cause discomfort.



- 2.2) Correctly fitted straps (as shown).  
To release straps, press the two prongs on the buckle sides.

## CLEANING & MAINTENANCE

### **Cleaning Instructions:**

Clean down with a damp cloth and a mild detergent.  
When wet allow to dry naturally in a warm dry atmosphere.  
**DO NOT** dry in any drying machine.  
**DO NOT** wash in an automatic washing machine.  
**DO NOT** iron.

### **Maintenance Instructions:**

**CHECK** for any signs of fraying and that buckles work correctly.  
**CHECK &** make sure that all fixings are secure & not damaged, especially after making any adjustments.  
**PERIODIC** cleaning will help prevent damage to buckle.