

# COSY TOES SELF ASSESSMENT FORM





So we can get your Cosy Toes exactly right please note the following information on the assessment form. If you need help or more information call us direct on (44) 01283 520400

Client Name: .....

Contact Phone Number: ..... Email: .....

Please state which unit of measurement you are using here: Please  box.  Inches  Millimetres



- A** Chest to front of knee  
(Measure in a straight line as shown by arrow A on the diagram) With backrest reclined: ..... Upright: .....  
(if applicable) (as pictured)
- B** Front of knee to toes  
(Measure in a straight line whilst feet are resting on the footrest / footbox (as shown by arrow B on the diagram) .....
- C** Front of toes to front of seating system  
(Measure in a straight line as shown whilst feet are resting on footrest / footbox (as shown by arrow C on the diagram) .....
- D** Seating system Overall width  
(Measure the outside edge of seating system to outside edge over the client's legs (as shown by arrow D on the diagram ) .....

### Armrests and Trays

If the client's chair has armrests please provide the measurement below. (As shown by arrow E on the diagram) The mounting points for a tray are usually attached to the armrests. If this is the case you may not be able to use Cosy Toes with a Wheelchair Tray. If you are unsure please contact us for advice.

- E** Armrest width  
(Outside edge to outside edge of the armrests (as shown by arrow E on the diagram) .....

**Please Turn Over for Colour Options**

If you are unsure or need advice, please contact us using the details below:

Tel: 0044 (0)1283 812860 Email: enquiries@specialorthotic.com Web: www.specialisedorthoticservices.co.uk

## Name options (free of charge)

Options: Please  box.

Embroidered Name: .....

## Lining Options

Please note there is an additional charge for Sheepskin.

Options: Please  box.

Fleece Liner       Sheepskin Liner

## Colour & Style

Options: Please  box.

Coloured Fabric\* With Black Trim

\*Select Colour & Material from options below:

## Pastel Polyester



Hot Pink



Purple



Lime



Crimson



Turquoise



Royal



Black



Orange

Please note that Cosy Toes is wipe clean only and is not machine washable.

If you are unsure or need advice, please contact us using the details below:

Tel: 0044 (0)1283 812860    Email: [enquiries@specialorthotic.com](mailto:enquiries@specialorthotic.com)    Web: [www.specialisedorthoticservices.co.uk](http://www.specialisedorthoticservices.co.uk)