



Equipoise Straddle Seat



Specialised
Orthotic
Services
By *drive*



Functional sitting
For therapy, work & play



Helps children with low levels of sitting ability work on sitting balance
from a young age.



specialisedorthoticservices.co.uk



Made in the UK

SOS08 ISS03 0518



Equipoise Straddle Seat

The Equipoise has been developed as an alternative to function seating. It helps children with low levels of sitting ability maintain a balanced seating position.

Conventional seating relies on a 90/90 sitting position which has basic intrinsic disadvantages whilst working on sitting balance. The Equipoise Straddle Seat uses the sitting position adopted when riding a horse. By reducing the hip flexion the client automatically produces a lordosis and puts the extensor muscle group at an advantage.

The straddle position minimises the effects of sitting with a poor centre of gravity and helps the child achieve a more balanced seating position.



Freestanding tray with base plinth



Base pads

Key features:

Available in 2 different versions

Equipoise is available with a mobile base or static nursery version. The static nursery version is available in 2 sizes.

Functional sitting for therapy, work & play

The straddle position helps the client sit up straight, enabling functional movement of the upper body. Ideal for short therapy sessions and longer work and play sessions.

Angle & height adjustable tray

Base pads to accommodate growing children

The nursery version of the Equipoise is available with detachable base pads. This allows the height to be adjusted quickly, ideal for groups of children.

Suitable for young children

Equipoise gives children with low levels of sitting ability the opportunity to work on sitting balance from a young age.

Removable covers for easy cleaning



1. Positioning Velcro Bands
2. Foot Support (height & angle adjustable)
3. Tray (angle adjustable)
4. Saddle Support
5. Mobile Tilt in Space Base
6. Static nursery version
7. Base pads



Specialised
Orthotic
Services
by *drive*

For more information, advice or to book a demonstration

call: 0044 (0)1283 520400

email us: enquiries@specialorthotic.com

visit our website: www.specialisedorthoticservices.co.uk

